### **Experiment 5**

## **EDTA Titration: Calcium in Calcium Supplements**

#### **Student Handout**

#### **Purpose**

To determine the amount of calcium in a calcium supplement tablet by EDTA titration.

#### **Background**

Calcium is an important element for our body. It is vital for the development of bones and teeth. More than 95% of calcium in our body can be found in bones and teeth. The long-term calcium deficiency leads to the decrease in bone density and strength. In fact, osteoporosis has affected many elderly people. To avoid this problem, calcium tablets may be taken to supplement the daily intake from diet.

The calcium supplements consist mainly of calcium salts. Other ingredients are binders, colouring agents, flavouring agents, etc. In the present experiment, the amount of calcium in a supplement tablet will be determined by titration with a standard solution of ethylenediaminetetraacetic acid (H<sub>4</sub>EDTA).

EDTA is commonly used as a chelating agent. It forms very strong 1:1 complex with most metal ions. The structure of  $H_4EDTA$  is shown below.

As both the calcium ions and the calcium-EDTA complex are colourless, a suitable indicator that will change colour at the equivalence point has to be employed. Calmagite is used in the present experiment. It binds less strongly to the calcium ions than EDTA and exhibits different colours in free and bind forms at certain pH value.

$$Ca^{2+}(In)(aq) + EDTA^{4-}(aq) \rightleftharpoons Ca(EDTA)^{2-}(aq) + In(aq)$$
 (wine red) (blue)

where *In* stands for indicator.

#### **Safety**

Handle all chemicals with great care. Avoid direct contact of chemicals with skin. Dispose of chemical waste, broken glassware and excess materials according to your teacher's instruction.

Safety information on the chemicals used in the investigation can be found in the Material Safety Data Sheet (MSDS). Consult your teacher for details.



### **Materials and Apparatus Available**

Calmagite indicator solution



Standard 0.05 M EDTA solution



3 M HCl solution





3 M NH<sub>3</sub> solution





NH<sub>4</sub>Cl/NH<sub>3</sub> buffer solution with Mg(EDTA)<sup>2</sup>





Burette Conical flasks Analytical balance Mortar and pestle Beakers
Volumetric flasks
Pipette
Measuring cylinders

#### **Experimental Procedure**

Photos of the experiment are available at http://www.chem.cuhk.edu.hk/ssc.htm.

Make sure you have recorded all the observations and data.

#### 1. Preparation of sample solution

Write down the information for the calcium supplement, especially the amount of calcium present in each tablet. Estimate the portion of tablet used for titration with ~25 cm<sup>3</sup> of 0.05 M EDTA solution. Grind the tablet into powder. Weigh accurately the amount of powder needed and put it into a 250-cm<sup>3</sup> conical flask. Add ~5 cm<sup>3</sup> of 3 M HCl solution in small portions slowly to dissolve the powder.

#### Titration with standard EDTA solution

Add 50 cm<sup>3</sup> of distilled water and 5 cm<sup>3</sup> of 3 M aqueous NH<sub>3</sub> solution to neutralise the final solution in Part (1). Add ~1 cm<sup>3</sup> of buffer solution to obtain a pH value of ~10. Add a few drops of Calmagite indicator solution. The solution should be in wine red colour. Titrate with standard EDTA solution with continuous swirling until the wine red colour turns into blue. Repeat the titration twice.

#### **Questions for Further Thought**

- EDTA is not a primary standard. Suggest a method to standardise the EDTA solution.
- The present technique is employed to determine the hardness of water. Briefly describe what the hardness of water is and why we need to determine it.
- EDTA is used in chelation therapy. Describe briefly this treatment and the principle behind.

#### References

- 1. A. D. Eaton, L. S. Clesceri and A. E. Greenberg, Standard Methods for the Examination of Water and Wastewater, 19th Ed., American Public Health Association, Washington, 1995, pp. 2-35 - 2-38.
- http://web.centre.edu/che/che131 lab/calcium.html

# Experiment 5 EDTA titration: calcium in calcium supplements

# **Laboratory Report Form**

Date:			
Title:			
Purpose:			
Data and Calculation:			
Duta and Carcumiton.			
Sample Information:			
Brand name: Amount of calcium per tablet (labele	d):	mg	
Weight of tablet: g			
Estimated weight of powder used:	g		
Concentration of EDTA solution:	M		
	Trial 1	Trial 2	Trial 3
Weight of powder (g)			
Final reading of burette (cm <sup>3</sup> )			
Initial reading of burette (cm <sup>3</sup> )			
Volume of EDTA solution used (cm <sup>3</sup> )			
Calcium per gram of powder			
Average (g)			
% deviation			
/o deviation		l	<u> </u>
Amount of calcium per tablet:	mg		

Conclus	sions:
<u>-</u>	
Questio	ns for Further Thought
	TA is not a primary standard. Suggest a method to standardise the EDTA ation.
Answer:	
	e present technique is employed to determine the hardness of water. Briefly cribe what the hardness of water is and why we need to determine it.
Answer:	

# Experiment 5 EDTA titration: calcium in calcium supplements

3.	EDT princ			chelati	on	therapy.	Describ	e	briefly	this	treatment	and	the
Ans	wer:												
	-												
	-												
	-												
	-												
	-												